

LIASYS 450 DURABLE CUVETTE SYSTEM



80 optical grade cuvettes in 4 integrated segments set on in a high capacity reaction wheel, are continuously tested for cleanliness to eliminate carryovers.

WBL tests (Water Blank Level) run by the on-board Quality Assurance Program monitor each cuvette and warn if any of the segments includes a cuvette that does not meet the minimum quality requirements; if needed, each cuvette segment can be easily replaced by the system operator.

The intuitive and simple graphical interface provides immediate system and operation update and feed-back; look durable cuvette status, alarms, available test results, sample already processed, reagent positions, reagent containers levels, time toward tests or work-list completion.



LIASYS 450 durable cuvette segment

DURABLE OPTICAL CUVETTES - QC MONITORED AND USER CHANGEABLE